

EZ-150CA Euthanex Isoflurane/CO₂ Combination System

Isoflurane/CO₂ combination system

Features:

- Designed for investigators that need to retrieve animals quickly
- Utilizing up to (4) induction chambers with the system
- Easily switch over from Isoflurane to CO₂ during the procedure
- Isoflurane controlled by a TEC 3 Vaporizer
- Built in CO₂ Flow meter

The EZ-150CA System is designed for small scale use by individual investigators in their laboratories. This system first delivers Isoflurane using a TEC 3 Vaporizer then can be switched over to deliver a controlled CO₂ gas flow rate. System can accommodate (1-4) induction chambers.

